



## Respiratory Screen

### Introduction

Most acute respiratory tract infections are caused by commonly encountered viruses. As respiratory viruses cause a number of overlapping signs and symptoms, it is difficult to determine the causative agent by clinical presentation alone. Traditionally, we have used multiplex Polymerase Chain Reaction (PCR) as a means of rapidly screening for respiratory pathogens, as screening is of great benefit both in guiding optimal patient treatment and in deciding whether to implement measures to control the spread of infection (*Sze Hwei Lee et al., 2019*).

We have now advanced our assay portfolio to account for an increased number of patient samples, and to ensure our assays are in-line with the latest nucleic acid technology. The NxTAG Luminex® Respiratory Pathogen Panel (RPP) assay allows us to offer a comprehensive respiratory screen for up to 19 viral and 3 bacterial targets (see below), which are vital in the diagnosis of respiratory illness. The RPP assay is performed via the Luminex® MAGPIX, a microfluidics-based molecular platform.

The RPP assay kit has undergone extensive in-house verification in accordance with our UKAS laboratory accreditation, and is suitable for a variety of sample types (see below). An outstanding feature is the minimal volume of sample required for a significant diagnostic yield (>200µl).

### VERIFIED SAMPLE TYPES\*

- Nasopharyngeal swabs
- Throat swabs
- Bronchoalveolar lavage (BALS)
- Nasal aspirates
- Tracheal aspirates
- Sputum

**\*Represents patient sample types verified in-house for use in the Luminex® RPP assay. Other sample types may be tested upon request.**

## VIRAL PATHOGEN TARGETS

- Adenovirus
- Influenza A H1, Influenza A H3
- Influenza B
- Parainfluenza type 1-4
- Respiratory syncytial virus A (RSVA),  
Respiratory syncytial virus B (RSVB)
- Rhinovirus
- Enterovirus
- Seasonal Coronaviruses:  
229E/OC43/NL63/HKU1
- Human Metapneumovirus (MPV)
- Human Bocavirus

## BACTERIAL PATHOGEN TARGETS

- *Legionella pneumophila*
- *Mycoplasma pneumoniae*
- *Chlamydia pneumoniae*

### Pathogens detected by the Luminex® Respiratory Pathogen Panel (RPP).

A broad coverage of respiratory viruses in a single test can be of great benefit, both in ensuring optimal patient treatment, and in deciding whether to implement measures to control the spread of infection. Early diagnosis can help clinicians avoid the use of inappropriate treatments that may be costly and/or have potentially harmful side effects. This is particularly valuable when dealing with immunocompromised patients, among whom many respiratory infections carry a high morbidity and mortality (*Andrea Cortegiani et al., 2018*).

When viruses are detected by the Luminex RPP assay, we are often able to provide further subtyping of the pathogen, such as adenovirus and enterovirus (subject to sufficient viral load). In addition, we offer diagnostic testing for a number of prokaryotic and eukaryotic organisms through our established RT-PCR assays, including *Streptococcus pneumoniae*, *Mycoplasma pneumoniae*, *Mycobacterium tuberculosis* and *Pneumocystis jiroveci*. Please visit our website for a complete list of all our assays.

The diagnostic power of the Luminex RPP screen is not only applicable to our diagnostic services, but can also be used for high-throughput processing of research samples (*Poornima Ramanan et al., 2017*). The sensitivity of the RPP screen allows us to detect multiple viruses per sample, which can be of particular interest when interpreting a patient's clinical presentation and exposure. We are able to generate large amounts of data for research projects both quickly and cost-effectively using this approach, which is particularly advantageous in studies involving large cohorts. Please contact us if you would like to discuss using the screen as part of your research project (email: [info@micropathology.com](mailto:info@micropathology.com)).

## Influenza A and B

Influenza A and B are classified within the *Orthomyxoviridae* family, and target the respiratory tract within human hosts. Upon infection, timely diagnosis of influenza viruses may be critical for immunocompromised patients, due to the severity of illness caused by these viruses. Unlike most respiratory viruses, effective vaccines and antiviral therapies targeting influenza are available. The latter are chiefly neuraminidase inhibitors which include oseltamivir (marketed under “Tamiflu”), that interfere with viral cell entry/release (*Xiao-Guang Li, 2022*). The Luminex® RPP assay is able to differentiate between influenza types A and B, and distinguish between A subtypes H1N1 and H3N2.

## Metapneumovirus (MPV) and Respiratory Syncytial Virus (RSV)

MPV and RSV fall within the *Pneumoviridae* family, and are primarily thought to cause relatively mild respiratory tract infections. However, in young, elderly, and immunocompromised populations they may cause severe symptoms, resulting in hospitalisation (*Leigh M Howard, 2021*). RSV is commonly associated with croup in infants, and is well recognised as a significant cause of childhood mortality (*Marco Del Riccio et al., 2023*). Moreover, the association between RSV and the development of pneumococcal pneumonia can result in poor respiratory sequelae, highlighting the importance of RSV diagnosis and treatment (*Tom Wilkinson et al., 2023*). Similarly, MPV is believed to causes symptoms comparable to those of RSV, and has been well-defined to cause ARDS within an ICU setting (*Loreto Vidaur, 2019*). It is therefore important to diagnose both RSV and MPV to mitigate further transmission. The Luminex® RPP assay is able to distinguish between RSV and MPV infection, and further differentiate between RSV type A (RSVA) and RSV type B (RSVB).

## Parainfluenza viruses (HPIV) 1-4

Parainfluenza viruses (1-4) are members of the *Paramyxovirus* family, and are responsible for a range of symptoms, including rhinorrhoea, tussis, bronchiolitis and pneumonia (*Akhil Chellapuri, 2022*). They are particularly important within young children, among whom they are one of the most common causes of croup and/or hospitalisation following RSV infection (*Ji Yoon Han, 2022*). All four parainfluenza viruses cause a full spectrum of respiratory symptoms, however croup is most commonly associated with type 1, while bronchiolitis and pneumonia are most often associated with types 1 and 2 (*Akhil Chellapuri et al., 2022*). The Luminex® RPP assay is able to distinguish between all four parainfluenza types (HPIV types 1-4).

## Rhinoviruses (HRV) & Enteroviruses (EV)

Both rhinoviruses and enteroviruses are part of the *Picornaviridae* family, *Enterovirus* genus. There are three observed rhinoviruses species, termed human rhinovirus (HRV) A, B and C, responsible for over half of all viral respiratory tract infections (*Aripuana Watanabe, 2010*). While rhinovirus infection is often mild, infection during early childhood has been extensively linked to asthma onset, as well as asthma exacerbation (*Yu-Tsun Su, 2020, Haiwen Liu, 2020*). There are four non-rhinovirus enterovirus family members, classified as EV-A, EV-B, EV-C & EV-D. Enteroviruses have a diverse tropism, and can cause a wide range of clinical syndromes including myocarditis, meningitis and paralysis. However, some enterovirus subtypes (notably EV-C and EV-D) are associated with respiratory system infection and can cause rhinovirus-like symptoms. Although typically mild, such symptoms can lead to wheezing, dyspnoea and pneumonia. In 2014, EV-D68 was identified as the cause of a notable outbreak of respiratory disease in children, linked to neurological tropism and hospitalisation. Currently, the incidence of EV-D68 infections continue to rise, attributed to both a mutation-driven increase in virulence, as well as advancements in molecular diagnostics and accurate identification (*Léna Royston, 2016*). Our in-house PCR assays allow us to distinguish between rhinovirus and enterovirus, and we are able to perform genotyping upon request.

## Coronaviruses

### Pandemic (SARS-Cov-2)

Micropathology Ltd offers an in-house assay for the detection of SARS CoV-2 RNA, which is available upon request, and is routinely performed when the viral respiratory screen is requested.

### Seasonal (HCoV-229E, HCoV-OC43, HCoV-NL63, HCoV-HKU1)

Seasonal coronaviruses belong to the *Coronaviridae* family, with notable species including HCoV-229E, OC43, NL63 and HKU1, and are responsible for the infection of the upper respiratory tract. Symptoms of seasonal coronavirus often include coryza, pyrexia and dyspnoea, and while most infection is self-limiting, hospitalisation may be necessary (*Taehee Kim, 2021, Zhi-Qi Zeng, 2017*). Due to the fastidious nature of coronaviruses, identification is primarily achieved using molecular methods such as the Luminex RPP screen, which can distinguish between HCoV-229E, OC43, NL63 and HKU1 with a high degree of both accuracy and sensitivity.

Please also note that this assay is not designed to detect Middle East respiratory syndrome coronavirus (MERS-CoV), due to the high level of containment required.

## Adenovirus (AdV)\*

Adenoviruses are classified within the *Adenoviridae* family, organised into 56 adenovirus serotypes and assigned to six species A to F, with certain serotypes associated with differing patterns of disease severity or clinical presentation. AdV is commonly associated with conjunctivitis, gastroenteritis and respiratory tract infections within humans, with respiratory manifestations including pharyngitis, rhinorrhoea, pyrexia, bronchitis and pneumonia (*Xinye Wang, 2021*). While AdV is one of several viral causes of acute respiratory disease syndrome (ARDS), infection may develop and impair alveolar gas exchange which, left untreated and undiagnosed, is frequently fatal (*Xuefei Chen, 2020*). Furthermore, Adenoviruses are of particular concern among hematopoietic stem cell transplant patients, in whom AdV may disseminate and increase both morbidity and mortality rates (*Marco W. Schilham, 2002*). Following Adenovirus detection, we are able to perform genotyping upon request.

## Bocavirus (HBoV)

Human Bocavirus is a member of the *Parvoviridae* family and was first identified in 2005. There are four recognised human Bocaviruses (HBoV 1-4), however HBoV1 has been associated with respiratory symptoms, and is frequently detected in the airways of children under 5yrs presenting with respiratory illnesses. However, understanding the pathogenicity of HBoV1 is complex as it typically co-associates with other respiratory viruses. Although there is a paucity of information, evidence suggests that HBoV1 has a causative role in acute wheezing and exacerbation of respiratory illnesses (*Yu Deng, 2012*). Rarely, life-threatening incidence of HBoV1 infection have been reported in infants (*Tina Uršič, 2012*). The Luminex® RPP assay is able to detect respiratory bocavirus (HBoV1-4).

## ***Legionella pneumophila (L. pneumophila, Lpn)***

*Legionella pneumophila* is a Gram-negative, facultative intracellular bacterium that targets and infects the lower respiratory tract, resulting in atypical pneumonia. While there are more than 58 distinct species of *Legionella spp.*, only 25 are known to be harmful to humans. The most prominent of the *Legionella* species is *L. pneumophila*, accounting for more than 95% of all *Legionella* cases (*James T. Walker, 2021*). The Luminex RPP panel detects only *Legionella pneumophila* within respiratory samples. However, since these organisms are notifiable to UKHSA, they will be appropriately caveated on the patient report as: 'Please note, *Legionella spp.* appears as a 'notifiable organism' causative agent on the UKGOV website and should be notified to relevant Public Health Authority'.

## ***Mycoplasma pneumoniae* (M. pneumoniae, Mpn)**

*Mycoplasma pneumoniae* is a Gram-negative bacterium, categorised into two major groups: P1 type 1 (P1-1) or P1 type 2 (P1-2), depending on the P1 protein genotype present. Infection with either type of *M. pneumoniae* is often mild and self-limiting, however patients can develop severe and fulminant disease without diagnosis and treatment (*Jasna Rodman Berlot, 2021*). Furthermore, since *M. pneumoniae* is a common cause of bacterial atypical pneumonia, rapid diagnosis and treatment will mitigate further clinical sequelae and transmission. The Luminex® RPP respiratory screen is able to detect *M. pneumoniae*.

Please also note that Micropathology Ltd offers an in-house PCR assay for the detection of *Mycoplasma pneumoniae*, which maintains a superior sensitivity to the respiratory screen.

## ***Chlamydia Pneumoniae* (C. pneumoniae, Cpn)**

*Chlamydia pneumoniae*, formally known as *Chlamydophila pneumoniae* (from 1999 to 2015), is a Gram-negative bacterium which targets and infects the upper respiratory tract within human hosts. *C. pneumoniae* primarily causes atypical pneumonia, with symptoms ranging from asymptomatic or mild to severe pneumonia, which can result in related complications such as asthma exacerbation, encephalitis and cardiovascular dysfunction (*Doriane Calmes, 2021, Hadeel A. Nasser, 2019*). For this reason, it is clinically important to pursue a diagnosis and treatment for *C. pneumoniae*. The Luminex® RPP respiratory screen is able to detect *Chlamydia pneumoniae*.

Please also note that Micropathology Ltd offers an in-house PCR assay for the detection of *Chlamydia pneumoniae*, which maintains a superior sensitivity to the respiratory screen.

## **Summary**

Our use of traditional PCR together with Luminex® provides an effective method of identifying the cause of respiratory tract infection and we are always happy to discuss tailoring the screen to suit clinical need. Assay performance is subject to regular internal quality testing and external quality assurance programmes in line with our UKAS laboratory accreditation. The respiratory screen has

proven to be highly useful as a research tool for high-throughput sample analysis and together with our genotyping assays can be used to provide a wealth of epidemiological data.

## References

Akhil Chellapuri, M. S.-W. et al, (2022). Human parainfluenza 2 & 4: Clinical and genetic epidemiology in the UK, 2013–2017, reveals distinct disease features and co-circulating genomic subtypes. *Influenza and Other Respiratory Viruses*, 16(6), 1122-1132. doi.org/10.1111/irv.13012

Andrea Cortegiani, F. M. et al, (2018). Immunocompromised patients with acute respiratory distress syndrome: secondary analysis of the LUNG SAFE database. *Critical Care*, 22(157). doi.org/10.1186/s13054-018-2079-9

Aripuana Watanabe, E. C. et al, (2010). Rhinovirus species and their clinical presentation among different risk groups of non-hospitalized patients. *Journal of Medical Virology*, 82(12), 2110-2115. doi.org/10.1002/jmv.21914

Doriane Calmes, P. H. et al, (2021). Chronic infection with Chlamydia pneumoniae in asthma: a type-2 low infection related phenotype. *Respiratory Research*, 22(72). doi.org/10.1186/s12931-021-01635-w

Hadeel A. Nasser, M. K. et al, (2019). Chlamydia pneumoniae Infection and Endothelial Dysfunction in Cardiovascular Disease. *Biochem. Cell. Arch.*, 19(1), 1057-1062. doi.10.35124/bca.2019.19.1.1057

Haiwen Liu, J. T. et al, (2020). Altered mast cell activity in response to rhinovirus infection provides novel insight into asthma. *Journal of Asthma*, 57(5), 459-467. doi.org/10.1080/02770903.2019.1585870

Jasna Rodman Berlot, U. K. et al, (2021). Mycoplasma pneumoniae P1 Genotype Indicates Severity of Lower Respiratory Tract Infections in Children. *Journal of Clinical Medicine*, e00220-21, doi.org/10.1128/JCM.00220-21

Ji Yoon Han, W. S. et al, (2022). Seasonal epidemiological and clinical characteristics of pediatric patients with human parainfluenza virus infection by serotype: a retrospective study. *Virology Journal*, 19(141). doi.org/10.1186/s12985-022-01875-2

Leigh M Howard, K. M. et al, (2021). Clinical Features of Human Metapneumovirus-Associated Community-acquired Pneumonia Hospitalizations. *Clinical Infectious Diseases*, 72(1), 108–117. doi.org/10.1093/cid/ciaa088

Léna Royston, C. T. et al, (2016). Rhinoviruses and Respiratory Enteroviruses: Not as Simple as ABC. *Viruses*, 8(1), 16. doi.org/10.3390/v8010016

Loreto Vidaur, I. T. et al, (2019). Human metapneumovirus as cause of severe community-acquired pneumonia in adults: insights from a ten-year molecular and epidemiological analysis. *Annals of Intensive Care*, 9(86). doi.10.1186/s13613-019-0559-y

Marco Del Riccio, P. S.-Y. et al, (2023). Defining the Burden of Disease of RSV in Europe: estimates of RSV-associated hospitalisations in children under 5 years of age. A systematic review and modelling study. *Biome*. doi.org/10.1101/2023.02.10.23285756

Marco W. Schilham, E. C. et al, (2002). High Levels of Adenovirus DNA in Serum Correlate with Fatal Outcome of Adenovirus Infection in Children after Allogeneic Stem-Cell Transplantation. *Clinical Infectious Diseases*, 35(5), 526–532. doi.org/10.1086/341770

Mila Prill, M. K. et al, (2012). Human Coronavirus in Young Children Hospitalized for Acute Respiratory Illness and Asymptomatic Controls. *The Pediatric Infectious Disease Journal*, 31(3), 235–240. doi:10.1097/INF.0b013e31823e07fe

Poornima Ramanan, A. L. et al, (2017). Syndromic Panel-Based Testing in Clinical Microbiology. *Clinical Microbiology Reviews*, 31. doi.org/10.1128/cmr.00024-17

Sze Hwei Lee, S.-Y. R.-C.-F.-Y.-R. et al, (2019). Performance of a multiplex PCR pneumonia panel for the identification of respiratory pathogens and the main determinants of resistance from the lower respiratory tract specimens of adult patients in intensive care units. *Journal of Microbiology, Immunology and Infection*, 52(6), 920–928. doi.org/10.1016/j.jmii.2019.10.009

Taehee Kim, H. C.-H. et al, (2021). Epidemiology and clinical features of common community human coronavirus disease. *Journal of Thoracic Disease*, 13(4), 2288–2299. doi:10.21037/jtd-20-3190

Tina Uršič, M. J. et al, (2012). Human bocavirus and other respiratory viral infections in a 2-year cohort of hospitalized children. *Journal of Medical Virology*, 84(1), 99–108. doi.org/10.1002/jmv.22217

Tom Wilkinson, S. B.-L. et al, (2023). Burden of respiratory syncytial virus in adults in the United Kingdom: A systematic literature review and gap analysis. *Influenza and Other Respiratory Viruses*, 17(9). doi.org/10.1111/irv.13188

Xiao-Guang Li, J. C.-J.-H.-Y.-B.-F.-Q. et al, (2022). Oseltamivir Treatment for Influenza During the Flu Season of 2018–2019: A Longitudinal Study. *Frontiers in Microbiology*, 13. doi.org/10.3389/fmicb.2022.865001

Xinye Wang, D. W. et al, (2021). Molecular typing of human adenoviruses among hospitalized patients with respiratory tract infections in a tertiary Hospital in Guangzhou, China between 2017 and 2019. *BMC Infectious Diseases*, 21, 748. doi.org/10.1186/s12879-021-06412-0

Yu Deng, X. G. et al (2012). High Viral Load of Human Bocavirus Correlates with Duration of Wheezing in Children with Severe Lower Respiratory Tract Infection. *Plos one*, 1(1). doi.org/10.1371/journal.pone.0034353

Yu-Tsun Su, Y.-T. L.-C.-S.-Y.-L.-I.-M.-C.-R.-C.-C. et al, (2020). High correlation between human rhinovirus type C and children with asthma exacerbations in Taiwan. *Journal of Microbiology, Immunology and Infection*, 53(4), 561–568. doi.org/10.1016/j.jmii.2018.12.001

Zhi-Qi Zeng, D.-H. C.-P.-Y.-X.-X.-S.-K. et al, (2017). Epidemiology and clinical characteristics of human coronaviruses OC43, 229E, NL63, and HKU1: a study of hospitalized children with acute respiratory tract infection in Guangzhou, China. *European Journal of Clinical Microbiology & Infectious Diseases*, 37(1), 363–369. doi:10.1007/s10096-017-3144-z



